



COMPATIBILITY

SPECIAL NEEDS PUSHCHAIRS & WHEELCHAIRS

Suitable for:

Make	Model	Size 51 or 70 brackets *	Adaptor Bar	Notes
Alvema	Ito	51		
Chunc	45	51		
Convoid	Cruiser	51	yes	
Fold a Way	Spring	51		
Hoggi	Bingo	51		
Invacare	Blade Plus	51		With small rear wheels & outrigged front only
Jazz EASyS	EASyS	70		
Lomax	Kid Active Folding & Recliner	51		
Mojo		51		
Ormesa	New Bug	51		
Patron	Corzo Xcountry	51		
	Corzino Xcountry	51		
	Corzino Basic	51		
	Tom 5 Streeter	70		
	Tom 4 X Country	51		
	Ben 4 X Country	51		
	Reha Triga	51		
	Tom 4 Xcountry CSI	51		
Ranger		51		
Rehabuggy, Otto Bock Modell	Kimba Spring Gr. 1 & 2	51		
Rehatec	Modell Mika	51		
Saturn	Transporter Buggy	51		
Swiftly	EASys	51		
Tendercare	Snappi	51		

Compromised fit

Maclaren	Major	51		Necessary location of mounting bracket prevents pushchair from being folded.
----------	-------	----	--	--

Not suitable for:

Activate	Tandem Buggy
Invacare	Action 3 Junior
Leckey	MyGo

Any wheelchair with large rear wheels that hide the pushchair/wheelchair chassis, so preventing attachment of mounting brackets

* Buggypod box includes one set of size 51 brackets. Size 70 brackets available separately from Buggypod stockists or direct from Revelo